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
CLINICAL DATA: Right breast lesions.

MACROSCOPIC:

Specimen received in formalin solution, in two containers, labeled as follow:

A) Right breast mass; consist of a piece of yellow fatty tissue, measuring: 3x2.5x1.5cm; that covered by skin measuring: 2.5x2x0.5cm; on external surface a lesion is seen, measuring: 1.5cm in diameter. Distance to margins is 0.5cm.

S.O.S: M/5 E: 90%

B) Right breast margins; re-excision, consist of a piece of yellow fatty tissue; measuring: 2.5x2.5x1.5cm; that covered by skin, measuring: 3x0.5x0.5cm.

S.O.S: M/4 E: 40%

Block description: 1) Medial margin 2) Superior margin 3) Lateral margin 4) Inferior margin

MICROSCOPIC:

A, B) See the diagnosis please.

DIAGNOSIS:

A) Right breast mass; excision:

- Ductal carcinoma in situ, high grade; solid, cribriform and comedo type; DIN-3.
- Inferior margin is closed to ductal carcinoma in situ, less than 0.1cm.
- Other surgical margins are free of tumor.
- Microcalcification are present in ductal carcinoma in situ.

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1: IDC Nuclear grade 2

- Ductal invasive neoplasia (DIN3) synonym of high-grade ductal carcinoma in situ
- Microcalcification are also seen.

2: Lymph node involved by carcinoma

3: Medial margin re-excision due to CDP diagnosis:

DCIS, high grade; cribriform and comedo type; DIN3 (pa-00-106)

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FROZEN SECTIONS REPORT:

- Involved medial and sub-areolar margins
- Free other surgical margins

IHC MARKERS (Block C):

Pan-CK: Positive in tumoral cells.


P63: Negative around some epithelial layers.

DIAGNOSIS:

A-D) Right breast mass and retroareolar, nipple and medial margins; partial mastectomy and margins re-excision:


- Invasive ductal carcinoma, NOS Type.
- Tumor size: 6cm in diameter. (pT3)
- Glandular differentiation: Score 3 (Less than 10% of tumor area is tubular formation)
- Nuclear pleomorphism: Score 2 (Large cell with vesicular nuclei, visible nucleoli and moderate variability in size and shape)
- Mitotic rates: Score 1 (About 3-4 mitoses per 10 hpf)
- Histologic grade: Grade 2 (According to Nottingham histologic score, total score 6/9)
- Sub areolar (Nipple) margin is involved by tumor.
- Other surgical margins are free of tumor after re-excision.
- Ductal carcinoma in situ, intermediate grade; solid, cribriform and comedo type; DIN-2; without extensive intraductal component.
- Microcalcification present in ductal carcinoma in situ and invasive component.
- Vascular and peri-neural invasion are noted.

c) A231-1: Medial Margin


A231-2: Medial Margin

A231-3: Inferior Margin

b) Subareolar (Nipple) margin detected by CDP

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CLINICAL DATA:
Received specimen in two separate formalin container labeled as:
1) Left partial mastectomy: consist of a piece of creamy yellow irregular soft tissue M: 8*6*5cm. On cut section an ill defined mass 0.8cm in diameter is seen. Distance of tumor to deep, inferior, lateral, medial superior are 1.5, 2.1, 0.1, 1.4 cm respectively. A small suspicious focus of tumor is seen in lateral margin.
SOS: 5/5 E: 5%
2) Left axillary lymph node, Sentinel node biopsy: Consist of a piece of creamy yellow irregular soft tissue TM: 6*4*3cm. SOS: 3/3 E: 60%

MICROSCOPIC:
See the diagnosis

DIAGNOSIS:
1) Left breast mass, Partial Mastectomy:
- Invasive ductal carcinoma
* nuclear grade: 2
* mitotic figure: 3
* Tubule formation: 3
* Size of tumor: 1cm
- The closest margin is lateral is <0.1cm
- All another margins are free at least 1cm
- Lymphovascular invasion not seen
- Perineural invasion are also seen
- Other area show fibrocystic changes

2) Left axillary lymph node, Sentinel node biopsy:
- One lymph node is involved by invasive carcinoma (1/2)

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A (2): IDC grade 2, columnar cell changes are also seen.
B (5): Fibrocystic change with moderate ductal hyperplasia (typical/atypical), sclerosing adenosis are also seen.
C (6): Non proliferating fibrocystic changes.
D (8): Fibrocystic change with moderate ductal hyperplasia, Apocrine metaplasia are also seen.
F (14): Unremarkable fatty tissue.

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- Free one Sentinel lymph node

DIAGNOSIS:
A, C and D) Right breast mass and superior margin and inferior margin; partial mastectomy and re-excision:
- Invasive ductal carcinoma, NOS Type.
- Tumor size: 6cm in diameter. (pT3)
- Glandular differentiation: Score 3 (less than 10% of tumor area is tubular formation)
- Nuclear pleomorphism: Score 3 (Large cell with vesicular nuclei, prominent nucleoli and marked variability in size and shape)
- Mitotic rates: Score 2 (About 6-7 mitoses per 10 hpf)
- Histologic grade: Grade 3 (According to Nottingham histologic score, total score 8/9)
- Superior, inferior and deep margins are involved in partial specimen and re-excised margins.
- Extensive stromal necrosis are seen about 50-70% of tumoral tissue.
- Microcalcification present in invasive component.
- Ductal carcinoma in situ not seen.
- Vascular invasion not present.
B) Sentinel lymph node; dissection:
- All five re-excision margins of tumor. (pN0)

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
a)

- Glandular differentiation: Score 2 (About 40% of tumor area is tubular formation)
- Nuclear pleomorphism: Score 2 (Large cell with vesicular nuclei, visible nucleoli and moderate variability in size and shape)
- Mitotic rates: Score 2 (About 6 mitoses per hpf)
- Histologic grade: Grade 2 (According to Nottingham histologic score, total score 6)
- All surgical margins are free of tumor after re-excision. (Infero margin is involved by tumor in frozen specimen)
- Ductal carcinoma in situ, high grade, comedo type: DIN-3; about 10-15% of tumoral area.
- Non tumoral tissue:
 - Extensive stromal fibrosis
 - Proliferative fibrocystic changes
 - Simple adenosis

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b)

- A (1): Fibrocystic change with columnar cell change.
- B (2): Non proliferating fibrocystic changes.
- C (4): Unremarkable fatty tissue.
- D (5): Invasive carcinoma with lobular feature.**
- F (8): Non proliferating fibrocystic changes.
- H (16): Non proliferating fibrocystic changes.
- I (17): Unremarkable fatty breast tissue.
- J (18): Unremarkable fibrofatty breast tissue.
- (19-22): fibro hemorrhagic fatty breast tissue and non proliferating - FCC with sclerosing adenosis. Suggestive for invasive carcinoma.**

After staining for Pan CK all epithelial cells within fatty tissue are stained so invasion is confirmed.

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
Clinical data
See our ref. No CM-99-2562

Specimen

- 1- Left breast areolar margin mass, lumpectomy, in fresh state for intra-operative pathology consult
- 2- Left axillary lymph nodes, excisional biopsy, in formalin
- 3- Right forearm mass, resection, in formalin

Macroscopy
Received specimens labeled as mentioned above in three (all samples and the specimens were transferred into 10% buffered formalin in lab, and the samples were fixed overnight with cold ischemia time of less than 1hr) and consist respectively as follows:

- 1- Part of breast tissue well-oriented by the surgeon with overlying skin (without nipple), without skeletal muscle strip attached, and the following gross description:
 - **Site and description of specimen:** Breast tissue measures 9x7x4 cm, superficial skin measures 3x1 cm, free of tumor grossly
 - **Tumor greatest diameter:** 2.7 cm
 - **Distance of tumor from surgical margins grossly:** Superior margin: 1cm; Inferior margin: 0.3cm; Medial margin: 0.2cm; Lateral margin: 3cm; Superficial margin & skin: 0.5cm; Deep margin: 2.5cm
 - **Other parts of breast tissue:** Grossly unremarkable
 - **Description of blocks:** Tumor: A1-3; Superior (inked black) and inferior (inked green) margins: A4; Medial (inked black) and lateral (inked green) margins: A5; Superficial (inked black) & deep (inked green) margins & skin: A6; FS (frozen section): Inferior margin fibrotic area

Frozen section report: The surgical margins are free

- 2- A fatty tissue piece measuring 5x3.5x2 cm, containing two lymph nodes, 1.5 & 2.5 cm in maximum diameter, embedded totally in three blocks, B1,2= One lymph node sliced, B3= Another lymph node sliced
- 3- A piece of homogeneous fatty tissue measuring 4.5x4x1.5 cm; SOS 5/3, E 30%; #C1-3

Microscopy
See the diagnosis please.

Immunohistochemistry
See the diagnosis below for the result of ER, PR, and Her2, IHC evaluation on tumor, performed on block #A3.

Dx:

- 1- Left breast areolar margin mass, lumpectomy (FS/PS):
 - **Invasive ductal carcinoma NST**

- Tumor size: 2.7 cm
- Histologic grade according to Nottingham system: Grade II/III (Tubule formation: 3, Nuclear grade: 3, Mitotic activity: 3; Total score: 9/9)
- Non-extensive high-grade DCIS present without comedo necrosis
- Lymph-vascular invasion is not identified in the sections examined
- Additional findings: Small fibroadenoma (at inferior margin)
- Stromal tumor-infiltrating lymphocytes: NI
- All surgical margins and skin are free of tumor
- Immunohistochemistry:
 - > HER2 (Clone SP3): Equivocal (2+): Incomplete and/or weak to moderate circumferential membrane staining within > 10% of the invasive tumor cells, or intense complete circumferential membrane staining \leq 10% of invasive tumor cells. HER2 FISH or CISH is recommended.
 - > ER (Clone SP1): Positive in 100% of invasive tumor cells with strong intensity
 - > PR (Novocasta PGR-312): Negative (Internal positive control present)

2- Left axillary lymph nodes, excisional biopsy:

- **Two found lymph nodes are free**

3- Right forearm mass, resection:

- Lipoma

ICD-O: C50.9, M8500/3

Pathologic staging (based on information available to the pathologist): pT2, pN0, pMX

Note: Report format based on protocol for examination, and template for reporting results of biomarker testing of specimens from patients with carcinoma of the breast (Feb 2020, College of American Pathologists).

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CLINICAL DATA: Known case of left breast carcinoma.

MACROSCOPIC:

Specimen received in fresh status and formalin solution, in six containers, labeled as follow:

A) Left breast mass; consist of a piece of yellow fatty tissue, measuring: 18x17x5cm; that covered by skin measuring: 10x6x0.5cm; on cut sections multiple lesions are seen; first lesion, measuring: 2cm in diameter; distance to anterior, inferior, deep and lateral margins are 3, 7, 3 and 16cm. Superior-anterior-medial margins of first lesion are involved by tumor. The other lesions, measuring: 0.5-2cm; distance to margins are 2-6cm.

S.O.S: M/15 E: 15%

Block description: 1) Lesion and superior-anterior-medial margins 2-4) Lesion 5-13) Random of in situ areas 14) Deep margin 15) Anterior margin and skin

B) Left Sentinel lymph node; consists of three lymph nodes, measuring: 1-2cm in diameter.

S.O.S: M/3 E: 100%

C) Left breast; superior margin, consist of two pieces of skin tissue; measuring: 8x5x2cm.

S.O.S: 2/2 E: 10%

D) Left breast; superior-medial margin, consist of a piece of yellow fatty tissue; measuring: 5x4x1.5cm.

S.O.S: 2/2 E: 10%

E) Left breast; inferior margin, consist of two pieces of yellow fatty tissue; measuring: 3x2x1cm.

S.O.S: 3/3 E: 60%

F) Left breast; medial margin, consist of a piece of yellow fatty tissue; measuring: 3x2x1cm.

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FROZEN SECTIONS REPORT:

Ductal carcinoma in situ

-Superior, anterior and medial margins are involved

-Free other surgical margins

Free one Sentinel lymph node

DIAGNOSIS:

A, C-F) Left breast mass and superior margin and superior-medial margin and inferior margin and medial margin; partial mastectomy and re-excision:

Invasive ductal carcinoma, NOS Type.

-Tumor size: multifocal, 0.2-0.4cm in greatest diameter. (pT1a)

Glandular differentiation: Score 2 (About 30-35% of tumor area is tubular formation)

Nuclear pleomorphism: Score 2 (Large cell with vesicular nuclei, visible nucleoli and moderate variability in size and shape)

Mitotic rates: Score 1 (About 4-5 mitoses per 10HPF)

Histologic grade: Grade 1 (According to Nottingham histologic score, total score 5/9)

All surgical margins are free of tumor after re-excision.

Skin is free of tumor.

Ductal carcinoma in situ, high grade, solid, cribriform and comedo type; DIN-3 with extensive intraductal component (More than 25%).

Negative for peri-neural and vascular invasion.

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No tumoral necrosis is identified.

B) Left Sentinel lymph node; resection:

All three Sentinel lymph nodes are free of tumor.

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LN35-1: Reactive lymph node hyperplasia (RLH).

LN35-2: Reactive lymph node hyperplasia (RLH).

LN35-3: Lymph node involved by very small focus of invasive carcinoma about 1mm (micro-invasion).

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a)

a)

Report of mammography on the last page
Dear Dr Saravani

Both breasts and axilla US & ELUS
Breast composition: Heterogeneous fibroglandular
Fibroglandular tissue of both breasts is thick and hyperechoic along with mottled appearance and prominent ducts indicative of benign breast changes.
No evidence of suspicious mass or architectural distortion is seen in the right breast.
Accessory nipple is seen in the right 5 o'clock near IMF.

Left breast:

- 8 o'clock subareolar and periareolar parts, compatible with palpable mass, a hypoechoic mass with irregular shape and margin containing foci of internal microcalcification with internal vascularity in color doppler exam and score 3 on elastography (27*15mm).
- 8-9 o'clock 3 cm from nipple, compatible with area of thickening sensation and microcalcifications on mammogram, a non-mass area of mottled hypoechogenicity (40*10mm, likely DCIS).

No evidence of suspicious stiffness or alteration in strain ratio is seen in elastography on the right breast.

Left axillary fossa: an (11*5mm) lymph node with diffuse cortical thickening: 3mm. Another lymph node (14*6mm) with lobulated border is also seen here (both: borderline, B4)

No pathologic LAP is seen in the right axillary fossa.

R/breast BIRADS II: benign findings.
L/breast BIRADS V: highly suggestive of malignancy, Biopsy was performed upon request.

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CLINICAL DATA: Known case of left breast carcinoma.

MACROSCOPIC:
Specimen received in fresh status, in two containers, labeled as follow:
A) Left breast mass; consist of a piece of yellow fatty tissue, measuring: 1.5x1.1x5cm; that covered by skin measuring: 10x2x0.5cm and nipple, measuring: 2x1.5x1cm; on cut sections a mass is seen at posterior part of nipple; measuring: 2.1cm in diameter. Distance to skin, deep, inferior, superior, lateral and medial margins are 1, 8, 1.5, 2.5, 1 and 12cm. Another lesion is also seen at medial part of mass, distance from margins are 1-8cm.
S.O.S: M/14 E: 5%

Block description: 1) Superior margin 2) Deep margin 3) Lateral margin 4) In situ lesion 5) Mass 6-11) Lesion area 12) Mass 13, 14) Nipple and skin


B) Left Sentinel lymph node; consists of three lymph nodes, measuring: 0.5-1cm in diameter.
S.O.S: M/4 E: 100%

MICROSCOPIC:
A, B) See the diagnosis please.


FROZEN SECTIONS REPORT:

- Invasive ductal carcinoma and DCIS.
- Free surgical margins
- Free three Sentinel lymph nodes

DIAGNOSIS:
A) **Left breast mass; partial mastectomy:**
- Invasive ductal carcinoma, NOS Type.

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Tumor size: 2.1cm in diameter. (pT2)
Glandular differentiation: Score 2 (About 45-50% of tumor area is tubular formation)
Nuclear pleomorphism: Score 2 (Large cell with vesicular nuclei, visible nucleoli and moderate variability in size and shape)
Mitotic rates: Score 2 (About 6-7 mitoses per 10HPF)
Histologic grade: Grade 2 (According to Nottingham histologic score, total score 6)
All surgical margins are free of tumor.
Skin and nipple are free of tumor.
Ductal carcinoma in situ, high grade, solid, micropapillary and comedo type; DIN-3.
In situ component size: 4cm in diameter.
No vascular and peri-neural invasion is identified.
Negative for tumoral necrosis.


B) Left Sentinel lymph node; dissection:
All three lymph nodes are free of tumor.

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

LN13-B1: Reactive lymph node hyperplasia (RLH).
LN13-B2: lymph node show very monotonous appearance.
LN13-B3: Reactive lymph node hyperplasia (RLH) with foci of hemorrhage and fat necrosis.
LN13-B4: Reactive lymph node hyperplasia (RLH).

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









-Free surgical margin, right breast.
-Closed deep margin in left breast.
-Free other surgical margins, left breast.

DIAGNOSIS:

A) Right breast mass; partial mastectomy:
- Invasive ductal carcinoma, NOS Type.
-Tumor size: 1cm in diameter. (pT1b)

-Glandular differentiation: Score 2 (About 25% of tumor area is tubular formation)
-Nuclear pleomorphism: Score 2 (Large cell with vesicular nuclei, visible nucleoli and moderate variability in size and shape)
-Mitotic rates: Score 1 (About 3 mitoses per 10 hpf)
-Histologic grade: Grade 1 (According to Nottingham histologic score, total score 5/9)
-All surgical margins are free of tumor.

-Ductal carcinoma in situ, intermediate grade, solid and cribriform type; DIN-2; without extensive intraductal component.
-Microcalcification present in ductal carcinoma in situ and invasive component.
-Negative for vascular and peri-neural invasion.

B, E) Left breast mass and deep margin, partial mastectomy and re-excision:
- Invasive ductal carcinoma, NOS Type.
-Tumor size: 2.2cm in diameter. (pT2)
-Glandular differentiation: Score 3 (less than 10% of tumor area is tubular formation)

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-Nuclear pleomorphism: Score 3 (Large cell with vesicular nuclei, prominent nucleoli and marked variability in size and shape)
-Mitotic rates: Score 2 (About 7-8 mitoses per 10 hpf)
-Histologic grade: Grade 3 (According to Nottingham histologic score, total score 8/9)
-All surgical margins are free of tumor after re-excision.
-Ductal carcinoma in situ, high grade, solid, cribriform and comedo type; DIN-3; without extensive intraductal component.
-Microcalcification present in ductal carcinoma in situ and invasive component.
-Vascular invasion seen.

C) Right axillary lymph node; resection:
-All three lymph nodes are free of tumor. (pN0)

D) Left axillary lymph node; dissection:
-All eight lymph nodes are involved by tumor. (pN2a)
-Maximum diameter of metastatic focus is 2.5cm.
-Extranodal extension present.

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Left breast:
LN91-3: Lymph node involved by metastatic carcinoma about 90% (D2.80%).
LN91-6: Lymph node involved by metastatic carcinoma about 50%.
LN91-7: Lymph node involved by metastatic carcinoma about 70%.
LN91-8: Lymph node involved by metastatic carcinoma about 30%.
LN91-9: Lymph node involved by metastatic carcinoma about 80%.
LN91-10: Lymph node involved by metastatic carcinoma about 30%.
LN91-11: Lymph node involved by metastatic carcinoma about 50%.

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
c)

Cancer
Diagnostic Probe
(CDP)

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conventional pathology: Towards preoperative and intraoperative lymph node staging,” *Int. J. Cancer*, vol. 136, no. 4, pp. 743–751, 2015, doi: 10.1002/ijc.28742.

[9] Z. S. Miripour et al., “Electrically guided interventional radiology, in-vivo electrochemical tracing of suspicious lesions to breast cancer prior to core needle biopsy,” *Biosens. Bioelectron.*, p. 112209, 2020.

[10] Z. S. Miripour et al., “Electrochemical tracing of Hypoxia glycolysis by carbon nanotube sensors, a new hallmark for intra-operative detection of suspicious margins to breast neoplasia,” ., p. e10236.

[11] R. Mahdavi et al., “Intraoperative pathologically-calibrated diagnosis of lymph nodes involved by breast cancer cells based on electrical impedance spectroscopy; a prospective diagnostic human model study,” *Int. J. Surg.*, vol. 96, no. October, p. 106166, 2021, doi: 10.1016/j.ijsu.2021.106166.

[1] J. M. Jorns *et al.*, “Intraoperative frozen section analysis of margins in breast conserving surgery significantly decreases reoperative rates: one-year experience at an ambulatory surgical center,” *Am. J. Clin. Pathol.*, vol. 138, no. 5, pp. 657–669, 2012.

[2] M. S. Sabel *et al.*, “Residual disease after re-excision lumpectomy for close margins,” *J. Surg. Oncol.*, vol. 99, no. 2, pp. 99–103, 2009.

[3] A. Kobbermann *et al.*, “Impact of routine cavity shave margins on breast cancer re-excision rates,” *Ann. Surg. Oncol.*, vol. 18, no. 5, pp. 1349–1355, 2011.

[4] M. Rahman and S. Mohammed, “Breast cancer metastasis and the lymphatic system (Review),” *Oncol. Lett.*, vol. 10, no. 3, pp. 1233–1239, 2015, doi: 10.3892/ol.2015.3486.

[5] J. Hotton et al., “Pre-operative axillary ultrasound with fine-needle aspiration cytology performance and predictive factors of false negatives in axillary lymph node involvement in early breast cancer,” *Breast Cancer Res. Treat.*, vol. 183, no. 3, pp. 639–647, 2020, doi: 10.1007/s10549-020-05830-z.

[6] J. G. Reiter et al., “Lymph node metastases develop through a wider evolutionary bottleneck than distant metastases,” *Nat. Genet.*, vol. 52, no. 7, pp. 692–700, 2020, doi: 10.1038/s41588-020-0633-2.

[7] M. D. Alvarado et al., “SentimagIC: A Non-inferiority Trial Comparing Superparamagnetic Iron Oxide Versus Technetium-99m and Blue Dye in the Detection of Axillary Sentinel Nodes in Patients with Early-Stage Breast Cancer,” *Ann. Surg. Oncol.*, vol. 26, no. 11, pp. 3510–3516, 2019, doi: 10.1245/s10434-019-07577-4.

[8] M. Winter, R. Gibson, A. Ruszkiewicz, S. K. Thompson, and B. Thierry, “Beyond

Cancer Diagnostic Probe (CDP)

Version V. 1.0
Nov 21, 2021

www.NanoElecHealth.com